

 **decon**  
*Creating Freedom*



 **E-MOVE**





## Now you no longer need to adapt to the products! E-Move is the market's first push rim-controlled drive unit that is adapted to your needs.

Now you no longer have to adapt to the products. E-Move offers a variety of options to expand the functionality of a manual wheelchair and is completely adaptable to your needs.

When using the pushrim, E-Move supports your starting movement in a powerful way. You save energy and easily overcome longer distances or roads and obstacles that you normally avoid. With the new cruise control function, you can go long distances without wearing your self out. E-Move creates personal freedom and mobility!

## Specifications

Engine	AC servo motor, 24V 110W x 2
Battery type	NiMH, 24V x 6,7Ah / Li-iON, 25V x 10,08Ah
Weight wheel	7 kg
Total weight	NiMh 17 kg / Li-ion 18 kg
Increased seat width	Max 1.8 cm/each side
Wheel dimensions	20", 22", 24", 25" eller 26"
Air pressure	6 Bar, 600 kPa
Speed	0–8,5 km/h
Max range (with fully charged battery)	Ca. 18 km with NiMH, Ca. 27 km med Li-iON
Fuse	30A (flat)
Max userweight	150 kg
Max gradient	6 grader
CE marking	In compliance with MDR 2017/745 class 1
Battery charger	100 to 240V AC, 50/60Hz
Power rating	37V 26A (while charging)
Charging time	Ca 3 h NiMH / ca 4 h Li-iON



**Pleasae, scan the QR code for more information**



## E-move


---

### Battery positions

- Wheel placement
- Under seat in bag
- Behind back in bag

### Lithium-ion battery

- Removable
- 2-year warranty
- IATA-approved



## Gyro function and speed control

---

E-Move can be programmed to suit your unique needs: Extra power (forward, backward), turning sensitivity, hemiplegia adaptation, push rim sensitivity, imbalance in arm strength, speed settings, and gyro function.

With a second driving mode (optional), you can easily switch between two different driving programs, ensuring the right support in any situation.

E-Move is compatible with almost all standard wheelchair models and is available in wheel sizes 20", 22", 24", 25", and 26".

E-Move is easy to remove, making it simple to transport in a car. The E-Move battery is also approved for air travel.



Please, scan the QR code for more information



Together we can create freedom!



Follow us:



**decon.se**

