







**Netti 4U, 4U CED, 4U CE Plus, 4U CEDS Young
MED3112**






9




10


1	MED2000	1	
2	100 400 008	2	
3	100 400 014	1	
4	405241	1	
5	XSM8X40	1	
6	XSM8ML	1	
7	XSM8X10	1	
8	XSM8BR	1	
9	MED3112-01	1	
10	MED3112-02	1	
11	HW084	1	
12	XSM18M	1	




4U CED
MED1431 + MED1439



1	MED1430	1	
2	MC6SM5X22	1	




1	XOC-48160-B0	1	
2	MED1507	1	
3	XSM5BR	1	
4	MC6SM5X14	1	



1	XSM6BRTG	1	
2	Xklamma2	1	
3	XSM6BR	1	
4	XSM6X20	1	




MED1252 + MED1643



1	XA2-2612K-00	1	
2	XA2-2612E-00	1	
3	XSM6BRTG	1	
4	MED1317	1	
5	XSM6BR	3	
6	MC6SM6X35	1	
7	XSM6X16	2	
8	MED1316	1	
9	XSM6ML	2	



M3	1,2 Nm
M4	3 Nm
M5	6 Nm
M6	10 Nm
M8	25 Nm
M10	50 Nm
M12	85 Nm



M16x1	30 Nm
M18x1	40 Nm